REF


ReOrder No.

550771	XCP-ORA® Arm & Ring, 1-pk
550772	XCP-ORA® Arm & Ring, 2-pk
550777	XCP-ORA® Arm & Ring, 6-pk
550773	XCP-ORA® Ring Only, 1-pk
550774	XCP-ORA® Arm only, 1-pk
550598	Endo Ring, 1-pk
550597	Endo Arm, 1-pk
550500	Rinn® Universal Digital Sensor Covers, 500-pk

## Rinn® XCP-ORA® One Ring & Arm Positioning System

Système de positionnement à un anneau et un bras / Sistema con un braccio e un anello di posizionamento / Ein-Ring- und Armpositionierungssystem / Sistema de posicionamiento con un anillo y brazo / Systeem met één positioneringsring en -arm / Jeden system ramienia pozycjonującego i pierścienia / Sistema de posicionamento de um anel e braço / **Ενιαίο σύστημα τοποθέτησης δακτυλίου και βραχίονα** / Yksi rengas- ja varsijärjestelmä / Èt ring- og armpositioneringssystem / Ett posijoneringsssystem for ring og arm / Ett ring- och armpositioneringssystem

For XCP®, XCP-DS®, XCP-QT® and XCP-DS FIT® Biteblocks

### Positioning Ring and Arm Kit

Kit de bras et anneau de positionnement / Kit braccio e anello di posizionamento / Positionierungsring- und arm-Kit / Juego de anillo y brazo posicionador / Kit met positioneringsring en -arm / Zestaw ramienia pozycjonującego i pierścienia / Kit de anel e braço de posicionamento / Kit τοποθέτησης δακτυλίου και βραχίονα / Kohdistusrengas ja -varsisarja / Ring- og armpositioneringssæt / Posisjoneringssett for ring og arm / Sats med positioneringsring och arm

Manufactured by:  
Dentsply Rinn  
1301 Smile Way  
York, PA 17404 USA  
www.dentsplysirona.com

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




### Alignment adjusted for biteblock flexing

Alignment ajusté pour la flexion du bloc de morsure / Allineamento regolato in funzione della flessione del biteblock / Ausrichtung aufgrund Bissblockbiegung angepasst / Alineación ajustada a la flexión del bloque de mordida / Uitlijning aangepast voor buiging van het biteblock / Dopasowanie dostosowane do odgięcia wafu zwarciowego / Alinhamento ajustado para flexibilidade do bloco de mordida / Ευθυγράμμιση προσαρμοσμένη για κάμψη του στοματοδιαστολέα / Kohdistus säädetty purentapalan venymistä varten / Vinkling tilpasset til bøining af bideblok / Innretning justert for bøying av bitepinne / Innretning justerad för böjning av bitblock

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

## ENGLISH

**Indications for use:** XCP® and XCP-ORA® arms and rings are indicated for positioning and aligning dental intraoral radiographs.

**Contraindications:** None known.

**Warnings:** Sterilize device before use. DO NOT use phenol-based glutaraldehyde or ultrasonic cleaners. DO NOT chemoclave, dry-heat sterilize or cold sterilize. Only sterilize with steam sterilization according to the provided sterilization instructions. DO NOT exceed 273°F (134°C). It is the user's responsibility to qualify alternative methods of disinfection and sterilization. Failure to follow properly validated sterilization processes and approved aseptic techniques may result in cross contamination. Any method of sterilization will shorten the life of plastic parts.

**Precautions:** Follow sterilization instructions carefully to prevent plastic damage during sterilization.

**Adverse Reactions:** None known.

**Residual Risks:** Failure to follow the validated reprocessing instructions provided may result in exposing the patient to cytotoxic residuals and/or cross-contamination.

### Step-by-Step Instructions:

1. Insert positioning arm into the appropriate color-coded channel of the aiming ring.
2. Attach appropriate biteblock for area to be imaged. Correctly assembled, the imaging area should be centered when viewed through the aiming ring.
3. Position holder in patient's mouth and slide ring close to patient's face.
4. Align x-ray tube to ring.
5. Take exposure.
6. Adjust as needed for the next desired position.

**Reprocessing:** Arm and Ring: Disassemble instrument(s) and separate all parts prior to reprocessing.

XCP-DS Fit™ Biteblocks: Separate the silicone band from the biteblock and retaining clip prongs, leaving the band loosely attached to the biteblock.

### Cleaning: Manual / Mechanical (Ultrasonic)

1. Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
2. Prepare an enzymatic instrument cleaning solution and immerse the device per the manufacturer's instructions.
3. Rinse the device under potable running water for at least 30 seconds and dry using a disposable lint-free cloth.
4. Visually inspect the device for visible soil. If soil is seen, repeat the cleaning procedure.

### Cleaning and Disinfection: Automated – Instrument Washer-Disinfector

Use only properly maintained, calibrated, and approved washer-disinfectors according to ISO 15883.

1. Rinse the disassembled device under potable running water while brushing with a soft-bristled brush to remove gross soil.
2. Follow manufacturer instructions for instrument washer solution, observing concentration rates and contact times as well as the instrument washer's recommended cleaning and/or thermal disinfection cycle (1 minute at 90°C).
3. Visually inspect the devices for visible soil and repeat the cleaning and disinfection procedure if necessary.

**Sterilization:** Ensure that the device is thoroughly dry prior to sterilization. Enclose the device into an appropriately sized sterilization pouch intended for steam sterilization. Place the pouched device into the sterilizer per the manufacturer's instructions. Sterilize the device according to one of the recommended methods listed below.

**Sterilize metal and plastic components in separate pouches to avoid damage to the components. DO NOT exceed 273°F (134°C).**

Steam Sterilization Cycle	Minimum Temperature	Minimum Full Cycle Exposure Time	Minimum Dry Time
Pre-vacuum	270°F (132°C)	3 min	16 min
Steam-Flush Pressure-Pulse	270°F (132°C)	3 min	16 min
Gravity Steam Sterilization	270°F (132°C)	10 min	15 min

**Alternate Method:** Place the non-pouched device into the sterilizer and sterilize according to one of methods listed above. Non-pouched device must be used immediately following sterilization.

**Storage:** To maintain sterility, the device should remain in the sterilization pouch until ready for use.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the Competent Authority of the Member State.

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© 2019 Dentsply Sirona