


endosuccess canal access prep

The micro-blades are less aggressive than diamond and their coating makes these tips very durable.

F88180



CAP1, CAP2, CAP3 tips, autoclavable metal support and universal wrench



Canal access preparation



CAP1

Micro-blade tip length 12mm, taper 6%

Active lateral part for:

- Finishing walls and polishing.
- Removing temporary cement and dentinal residues.
- Removing dentin overhangs.

Non-active end to prevent the risk of perforating the pulp chamber floor.



CAP2

Micro-blade tip, length 9mm, taper 5%

Active lateral part and extremity used by sweeping method to remove dentine bridges.

- Location of the MB2 (2nd mesiobuccal canal) and search for hidden canals.
- Preparation of the pulp chamber.
- Removal of the dentine layer which may hide the access to the MB2 canal.



CAP3

Micro-blade tip, length 8mm, taper 6%

The CAP3 tip has a very pointed extremity indicated for:

- Locating and opening the calcified canals.
- Fragmenting calcifications or pulp stones in the pulp chamber.
- Loosening fiber posts.
- Locating accessory canals.

Due to its very sharp point, the CAP3 tip must be handled with care (visual aids recommended).



ET18D

Diamond-coated steel tip 76µm, length 18mm, taper 5%

- Finishing the access cavity.
- Removing dentine overhangs, calcifications and filling materials.



ETBD

Diamond-coated ball tip, length 20mm, taper 5%

Searching for canals and locating calcified canals.

irrisafe



IRRI20,25

Irrigation



IRRISAFE

Passive ultrasonic irrigation (PUI) files of different lengths and diameters

Irrisafe™ safely* eliminates the smear layer, dentine debris and bacteria from the root canal. Its blunt tip prevents any risk of perforating the apex or the canal walls.

Irrigation once the root canal has been prepared.

- 20ml of irrigant (NaOCl) are injected into the canal.
- Irrisafe™ is inserted 2mm short of the working length and activated by performing withdrawal movements to flush the debris and the smear layer upwards.
- Repeated 3x 1 minute in each canal.



K FILES

Files of different lengths and diameters, taper 2%

Irrigation, withdrawal of calcified dentine and gutta percha, and withdrawal of broken instruments.

For irrigation ultrasonic files are used with a disinfectant solution. To provide a final decontamination, use sodium hypochlorite until the smear layer is removed.

K files are very sharp instruments and should be handled with precision. However they are flexible and can therefore be pre-bent.

K10, 15, 25, 30
FILES

* Van Der Sluis L.W.M. Passive ultrasonic irrigation of the root canal: a review of the literature. Int. Endodont. J. 2007; 40; 4: 415-428

CAP1
F88181CAP2
F88182CAP3
F88183ET18D
F88017ETBD
F88020