## ENGLISH

### Vinylpolysiloxane Impression Material System

- Type 3: low consistency - light bodied
- Type 2: medium consistency - medium bodied
- Type 1: high consistency - heavy bodied
- Type 0: very high consistency - putty

### General Description

Extrude is a system of vinylpolysiloxane impression material viscosities suitable for all crown and bridge, edentulous and implant impressions. Extrude is offered in four viscosities to address the complexity of each case.

### Inserting Cartridge into Extruder Gun

- (Note: Extruder gun comes fully assembled and requires no assembly by the user)
- Push metal lever at back of extruder up (toward plunger assembly) and pull back on plunger assembly until assembly is fully retracted.
  - Lift up on lock assembly on top of extruder exposing slot for insertion of cartridge.
  - Insert cartridge with notch in collar of cartridge down to match protrusion of extruder.
  - Snap lock assembly down locking cartridge into place.

### Attaching Mixing Tip to Cartridge

- While holding cartridge with cap toward you, rotate cap 90° counter clockwise, moving alignment tab on cap to align with notch on cartridge.
- Separate cap from cartridge by pulling cap off.
- For first use bleed cartridge by dispensing a small amount of material to ensure even extrusion from each side of the cartridge. Wipe excess material from cartridge.
- Take the mixing tip and align alignment tab with notch on cartridge. Push mixing tip down onto cartridge. Rotate 90° clockwise to lock mixing tip onto cartridge.

- Attaching Mixing Tip to Mixing Tip**
- Insert the round end of the intra-oral tip into mixing tip.
  - Press intra-oral tip into mixing tip. A click should be heard when the intra-oral tip locks into place.

### Extruder Operation

- Squeeze Extruder trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released.
- For subsequent usage, remove the used mixing tip and squeeze the Extruder trigger to extrude material from the two openings of cartridge. Wipe end of cartridge. Install a new mixing tip and twist one quarter turn to lock in position.

- To remove the cartridge, hold the release lever up and withdraw plunger in a backward direction.

### Putty Instructions

- Brush a thin coat of Kerr VPS tray adhesive on the selected tray. Allow a minimum of 10 minutes to air dry for maximum adhesion.
- Measure equal amounts of Extrude XP Putty base and catalyst.
- Mix base and catalyst with finger tips until a homogenous color is achieved (approx. 30 seconds).
- Place the mixed putty into the adhesive coated tray.

- EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.**

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#### Model Pouring

- The impression may be poured 20 minutes after start of mix if a vacuum mixed die stone is used. If hand mixed die stone is used, delay pouring of impression for 1 hour. Pouring of impression may be delayed up to 14 days without loss of accuracy.

#### Storage

- Store at ambient temperature.

#### Shelf Life

- Refer to outer package for expiration date. Do not use after expiration date.

#### Disinfection of Extrude Impressions

- Extrude impressions may be safely disinfected in 2% glutaraldehyde solutions. Refer to manufacturer's directions for time, dilution and solution temperature recommendations.

#### Electro Plating of Extrude Impressions

- Extrude impression may be electro plated. Recommended plating solutions are as follows:

### Attaching Intra-Oral Tip to Mixing Tip

- Insert the round end of the intra-oral tip into mixing tip.
- Press intra-oral tip into mixing tip. A click should be heard when the intra-oral tip locks into place.

### Extruder Operation

- Squeeze Extruder trigger with moderate and even pressure to begin mixing material. Material will stop flowing after trigger is released.
- For subsequent usage, remove the used mixing tip and squeeze the Extruder trigger to extrude material from the two openings of cartridge. Wipe end of cartridge. Install a new mixing tip and twist one quarter turn to lock in position.
- To remove the cartridge, hold the release lever up and withdraw plunger in a backward direction.

### Putty Instructions

- Brush a thin coat of Kerr VPS tray adhesive on the selected tray. Allow a minimum of 10 minutes to air dry for maximum adhesion.
- Measure equal amounts of Extrude XP Putty base and catalyst.
- Mix base and catalyst with finger tips until a homogenous color is achieved (approx. 30 seconds).
- Place the mixed putty into the adhesive coated tray.

- EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.**

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