Dental Bone & Hard Tissue Surgical Instruments

Bone Scraper, Dental Bone Curettes, Bone Ronguers, Periotome

en-English

Attention

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

Intended Use:

These instruments are intended for various oral and maxillofacial surgical procedures, including bone preparation, contouring, debridement, and extraction assistance. They facilitate precise tissue manipulation while minimizing trauma and preserving surrounding anatomical structures.

The Instrument Specific Intended Use along with their variants is given below:

Bone Scrapers:

These are intended to be used for the harvesting and collection of autogenous bone material for use in bone grafting procedures. They assist in the controlled removal of cortical bone while preserving bone vitality to support bone augmentation.

Dental Bone Curettes:

These are intended to be used for the removal of granulation tissue, pathological soft tissue, and debris from bone surfaces and extraction sites during oral and periodontal surgical procedures. They assist in bone debridement and preparation for grafting or implant placement.

Bone Rongeurs:

These are intended to be used for cutting, trimming, and removing bone tissue during oral and maxillofacial surgical procedures. They assist in bone contouring, reshaping, and the

removal of small bone fragments to prepare surgical sites for implants, grafts, or extractions.

Periotomes:

These are intended to be used for severing the periodontal ligament to facilitate atraumatic tooth extraction while preserving surrounding bone structure. They assist in minimizing bone loss and preparing the socket for implant placement or prosthetic restoration

 Our Dental Bone & Hard Tissue Surgical Instruments are intended to be used by qualified dental professionals only in a professional dental facility.

Patient population:

These instruments are intended for use on adult patients undergoing oral and maxillofacial surgical procedures, including bone preparation, contouring, debridement, and extraction assistance. The appropriate instrument type and size must be selected based on the patient's bone condition, anatomical structure, and procedural requirements.

For Available Configurations/Models please visit the Oral & Microsurgery section in our Catalogue

Features

Bone Scrapers:

- Specifically designed for collecting high-quality autogenous bone material with osteogenic potential, supporting procedures like ridge augmentation and sinus lifts.
- Sharp, serrated blade allows precise bone removal with minimal trauma to surrounding tissues, reducing patient discomfort and promoting faster healing.
- Ergonomic handle and blade design enable controlled scraping, ensuring optimal bone particle size and volume for grafting while minimizing wastage.
- Manufactured from surgical-grade stainless steel, ensuring sterility and durability for repeated use after proper sterilization.
- Versatile instruments suitable for implantology, periodontal regeneration, and maxillofacial reconstruction.

Dental Bone Curettes:

- Small, spoon-shaped or looped working ends with sharp, angled tips allow precise removal of granulation tissue, debris, or diseased bone from surgical sites.

- Made from high-grade stainless steel or titanium, ensuring durability, corrosion resistance, and ease of sterilization.
- Available in multiple sizes to accommodate both delicate and broad surgical areas, enhancing procedural accuracy and efficiency.
- Dual-ended design reduces the need for multiple instruments, improving workflow and minimizing procedural time.
- Sharp edges and precise contours facilitate thorough cleaning of bony defects, promoting better healing and reducing infection risks.

Bone Rongeurs:

- Specifically designed for cutting, trimming, and removing small pieces of bone during procedures such as alveoloplasty and extraction site preparation.
- Sharp, curved, or angled jaws provide controlled access to hard-to-reach areas, enabling precise bone removal.
- Some variants feature a spring-loaded hinge mechanism, allowing smooth operation and reducing hand fatigue during prolonged use.
- Manufactured from high-quality stainless steel, ensuring corrosion resistance, durability, and compatibility with sterilization processes.
- Available in different sizes and jaw configurations (straight, curved, or double action) to accommodate various surgical needs.

Periotome:

- Periotomes are designed with thin, sharp blades that allow precise severing of periodontal ligament fibers, minimizing trauma to the surrounding bone and soft tissues during tooth extraction or implant preparation.
- These instruments are specifically engineered to apply controlled force, reducing the risk of alveolar bone fractures compared to traditional extraction techniques.
- These instruments are less effective in cases of dense or highly calcified bone, as their thin blades may not provide sufficient leverage or strength for such conditions.
- They are ideal for atraumatic extractions, particularly in cases where preservation of the alveolar bone is critical, such as for immediate implant placement.

Operating instructions

   These instruments are non-sterile and must be cleaned and sterilized before first use.

   Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

Bone Scraper:

- Hold the instrument at a stable angle against the cortical bone surface.
- Apply controlled pressure while scraping in a back-and-forth motion to collect bone particles.
- Avoid excessive force to prevent unnecessary damage to surrounding tissues.
- Transfer collected bone material to the grafting site or storage container immediately after collection.

Dental Bone Curettes:

- Select an appropriately sized curette for the surgical site.
- Position the curette against the bone surface and use controlled, firm strokes to remove granulation tissue and debris.
- Maintain proper angulation for optimal tissue debridement without causing trauma.
- Regularly clean the curette tip during the procedure to ensure efficiency.

Bone Rongeurs

- Hold the instrument firmly with a stable grip.
- Position the jaws around the target bone area and apply controlled pressure to trim or remove bone fragments.
- Use a gradual closing motion to ensure precise cutting without excessive force.
- Avoid overloading the instrument to prevent damage to the cutting edges or handle mechanism.

Periotome:

- Position the periotome blade against the periodontal ligament space.
- Apply gentle, apical pressure while working around the tooth circumference.
- Maintain a stable grip and avoid excessive lateral force that may bend the blade.
- Continue ligament separation until the tooth is sufficiently mobilized for extraction.

Adverse Event

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

Cleaning & Sterilization

Cleaning

 Clean the device with care after each application to avoid the drying of contaminants.

Initial Cleaning

1. Remove any visible debris from the instrument using moisten a cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Use a soft brush and enzymatic cleaner to remove contaminants thoroughly.

Manual Cleaning

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

Automated Cleaning

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.
4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERILIZATION PROCEDURE.

 Ultrasonic cleaning may be performed in accordance with hospital or facility protocols. Ensure proper rinsing and drying after the process to prevent residue buildup.

Sterilization

Note: Ensure that the instruments are completely dry before sterilization.

Steam Autoclaving

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. In case of bone ronguers ensure that the hinges are lubricated with an appropriate medical-grade lubricant after sterilization to maintain smooth operation
5. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
6. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.
7. For bone ronguers, periodically check the hinge movement to ensure continued smooth function.

Important Notes

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.
- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturers, are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.
- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.

Disposal

Our Bone Scraper, Dental Bone Curettes, Bone Ronguers and Periotomes contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before

disposition. Please follow the local and national regulations or healthcare facilities' defined disposal and waste management policies.

Disclaimer

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

Explanation of utilized symbols

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)
	Manufacturing Date
	Manufacturer
	Authorized representative in the European Community
	Importer
	Keep in a cool, dry place
	Caution, read instructions

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