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tr

Instructions for Use

PROLENE™ Polypropylene (Monofilament) Suture

Sterile Synthetic Non-Absorbable Surgical Suture

Description

PROLENE™ Polypropylene Suture is a sterile, monofilament, synthetic, non-absorbable, surgical suture composed of ≥99.45 weight percent of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin.

The empirical molecular formula of the polymer is $(C_3H_6)_n$. PROLENE™ Suture is available undyed and pigmented blue with up to 0.55 weight percent phthalocyanine blue (Color Index 74160) to enhance visibility in the surgical field.

PROLENE™ Suture is available in a range of gauge sizes and lengths attached to needles of various types and sizes, and in presentations as described in the How Supplied section.

PROLENE™ Suture complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Non-Absorbable Monofilament Sutures and United States Pharmacopeia (USP) for Non-Absorbable Surgical Suture, except for a slight oversize in gauge size 7-0 (metric 0.5) (maximum oversize 0.007 mm). The European Pharmacopoeia recognizes units of measure Metric and Ph. Eur. sizes as equivalent that is reflected on the labeling.

PROLENE™ Suture is intended for use only by healthcare professionals who are trained in surgical suturing techniques.

Clinical benefits expected from soft tissue approximation and/or ligation are to promote wound healing by primary intention and to avoid bleeding from vessels during and after surgical procedure. A summary of safety and clinical performance can be found at the following link (upon activation):

<https://ec.europa.eu/tools/eudamed>.

Indications / Intended Use

PROLENE™ Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurosurgical procedures.

Contraindications

None known.

Warnings

Healthcare professionals should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing PROLENE™ Suture for wound closure, as risk of tissue separation / wound dehiscence leading to impaired healing may vary with the site of application and the suture material used.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Like all foreign bodies, PROLENE™ Suture may potentiate infection.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination that may lead to infection or transmission of blood-borne pathogens to patients and healthcare professionals.

Discard unintentionally opened / partially used / used devices and packages.

Precautions

As with any suture material, adequate knot security requires the standard surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the healthcare professional. The use of additional throws is particularly appropriate when knotting monofilament sutures. Improper knotting or damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

In handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.

Healthcare professionals should exercise caution when handling surgical needles to avoid inadvertent needle-stick injury that may result in transmission of blood-borne pathogens from contaminated needles. Broken needles may result in extended or additional surgeries or residual foreign bodies. In the event of a product malfunction before use, such as a bent, broken or detached needle, or suture damage, the product should be discarded and a new one obtained to begin the procedure. In the event of product malfunction during use, it is up to the discretion of the healthcare professional whether to continue or discontinue usage of the product and how to complete the procedure.

When attached to an EVERPOINT™ Cardiovascular Needle, note that EVERPOINT™ Needles are non-magnetic, and therefore are not compatible with magnetic devices such as counting trays.

Adverse Reactions / Undesirable Side Effects

Adverse reactions associated with the use of this device include tissue separation / wound dehiscence leading to impaired healing. Other associated adverse events include calculus formation in urinary or biliary tracts when prolonged contact with salt solutions such as urine or bile occurs, minimal inflammatory tissue / foreign body reaction, and transient local irritation at the wound site. Like all foreign bodies, PROLENE™ Suture may potentiate infection.

Suture breakage may result in hemorrhage. Suture material cutting or tearing through the tissue at the time of suture tensioning may cause soft tissue injury. Improper knotting or

damage to suture during use can cause extended surgery time or treatment failure and additional surgery.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

Healthcare professionals should convey adverse reactions, undesirable side effects, and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal postoperative course.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI) / Carcinogenic, Mutagenic, Toxic to Reproduction (CMR) / Endocrine Disrupting (ED) Safety Information

PROLENE™ Suture is MR safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360), based on human evidence and animal studies.

Application / Instructions for Use

Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique, and wound characteristics.

Discard needles in “sharps” containers. Discard unintentionally opened / partially used / used devices and packages according to your facility’s policies and procedures concerning biohazardous materials and waste.

Performance / Actions

PROLENE™ Suture elicits an initial, minimal inflammatory reaction in tissues that is followed by a gradual encapsulation of the suture by fibrous connective tissue. Implantation studies in animals show no significant change in the retention of tensile strength of the suture occurs during the lifetime of the implantation.

PROLENE™ Suture is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. Due to its relative biological inertness, it is recommended for use where the least possible suture reaction is desired. PROLENE™ Suture remains for the lifetime of the patient unless removed.

Sterility

PROLENE™ Suture is sterilized by ethylene oxide. Do not resterilize. Do not use if package is opened or damaged.

Storage

No special storage conditions required. Do not use after expiry date.

How Supplied

Please note that not all sizes are available in all markets. Please contact your local sales representative for size availability.

PROLENE™ Suture is available sterile, in sizes 10-0 through 2 (metric sizes 0.2 - 5.0) in a variety of lengths, with and without permanently attached needles.

PROLENE™ Suture is also available attached to CONTROL RELEASE™ Needles that enable the needles to be pulled off instead of being cut off.

PROLENE™ Suture is also available in configurations containing the following:


1. As HEMO-SEAL™ Needle Sutures combination in which the diameter of the suture and the needle wire have been more closely aligned to reduce the degree of needle-hole bleeding.
2. Pledgets, composed of 100 percent PTFE (polytetrafluoroethylene) felt, that are used as a buttress under sutures when tearing of friable tissue is a possibility.

PROLENE™ Suture is available in one, two or three dozen units per box.

Traceability

The following specific information can be found on the device packaging label: Catalogue number, Batch Code, expiry and manufacturing date, manufacturer name, address and the website and a Unique Device Identification bar code with the Global Trade Item Number information.

Symbols Used on Labeling

	Catalogue number		Undyed - Non-Absorbable - Monofilament - Uncoated - Suture
	Medical Device		Dyed - Non-Absorbable - Monofilament - Uncoated - Suture
	Do not use if package is damaged and consult instructions for use		Firm square cornered pledget threaded on suture
	Caution		Soft square cornered pledget threaded on suture
	Do not re-use		Firm round cornered pledget threaded on suture
	Do not re-sterilize		CONTROL RELEASE™
	Single sterile barrier system with protective packaging inside		Multi-strand
	Sterilized using ethylene oxide		
	Batch code		
	Date of manufacture		
	Use-by date		
	Unique Device Identifier		
	Manufacturer		
	Authorized representative in the European Community/European Union		
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.		
	Packaging unit		
	TEAR HERE		
	<p>Consult instructions for use or consult electronic instructions for use.</p> <p>EU: Call paper on demand help desk to get paper copies free of charge within 7 days.</p>		
	MR Safe		

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AB: Basılı kopyaları 7 gün içinde ücretsiz olarak talep üzerine almak için danışma hattını arayın.
Ethicon, LLC
Ethicon LLC
Johnson & Johnson Medical GmbH