elite HD+




A-Silicone dental impression material

putty soft
tray material
monophase
regular body
light body
super light body

IT - EN - FR - ES - DE - PT - NL - FI - EL - DA - SV
NO - RU - CS - RO - HR - HU - SK - LV - LT - ET - BG


IMPRESSION TAKING PROCEDURE
Two-phase technique (double impression):

1.1

1.2

1.3



1.4


1.5

One-step technique (simultaneous):

2.1




2.2





2.3

PREPARING THE DEVICES: STEP-BY-STEP INSTRUCTIONS
250 / 450 ml jars (Elite HD+ Putty Soft)







90 ml tubes (Elite HD+ Regular Body / Elite HD+ Light Body)








1:1 50 ml cartridges


**(Elite HD+ Tray Material / Elite HD+ Monophase / Elite HD+ Regular Body
 Elite HD+ Light Body / Elite HD+ Super Light Body)**

5:1 380 ml cartridges (Elite HD+ Maxi Putty Soft / Elite HD+ Maxi Tray Material)

**DISINFECTING AND CLEANING THE DEVICE:
GENERAL STEP-BY-STEP INSTRUCTIONS**

ATTENTION: For correct manual cleaning and disinfection of the cartridge and mixing tip read sections 12.5.2 and 12.5.3 carefully

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A-silicone dental impression material.

2. PRODUCT DESCRIPTION

Elite HD+ Putty Soft (mixing ratio 1:1 and 5:1): **high viscosity** vinyl polysiloxane, recommended for the two-phase (double impression) technique in combination with Elite HD+ Light Body or Elite HD+ Super Light Body fluids, or the one-step (simultaneous) technique in combination with Elite HD+ Regular Body or Elite HD+ Light Body fluids.

Elite HD+ Tray Material (mixing ratio 1:1 and 5:1): **high viscosity** vinyl polysiloxane, recommended for the two-phase (double impression) technique in combination with Elite HD+ Light Body or Elite HD+ Super Light Body fluids, or the one-step (simultaneous) technique in combination with Elite HD+ Regular Body or Elite HD+ Light Body fluids.

Elite HD+ Monophase (mixing ratio 1:1): **medium viscosity** vinyl polysiloxane, recommended for the monophase technique as a single viscosity material.

Elite HD+ Regular Body (mixing ratio 1:1): **medium viscosity** vinyl polysiloxane, recommended for the one-step (simultaneous) technique, in combination with Elite HD+ Putty Soft or Elite HD+ Tray Material.

Elite HD+ Light Body (mixing ratio 1:1): **low viscosity** vinyl polysiloxane, recommended for the two-phase (double impression) technique or the one-step (simultaneous) technique, in combination with Elite HD+ Putty Soft or Elite HD+ Tray Material.

Elite HD+ Super Light Body (mixing ratio 1:1): **low viscosity** vinyl polysiloxane, recommended for the two-phase (double impression) technique, in combination with Elite HD+ Putty Soft or Elite HD+ Tray Material.

3. PACKAGING

- 250 base jars + 250 ml catalyst (Elite HD+ Putty Soft)

- 450 base jars + 450 ml catalyst (Elite HD+ Putty Soft)

- 380 ml cartridge (Elite HD+ Maxi Putty Soft; Elite HD+ Maxi Tray Material)

- 50 ml cartridge (Elite HD+ Tray Material; Elite HD+ Monophase; Elite HD+ Regular Body; Elite HD+ - Light Body; Elite HD+ Super Light Body)

- 90 ml tube (Elite HD+ Regular Body; Elite HD+ Light Body)

4. COMPOSITION

Elite HD+ Putty Soft Normal / Elite HD+ Maxi Putty Normal / Elite HD+ Putty Soft Fast / Elite HD+ Maxi Putty Fast: vinyl polysiloxane, polydimethylsiloxane, silica fillers, methylhydrogenosiloxane, hydrocarbons, pigments, inorganic aluminium fillers, additives, organoplatinum complex.

Elite HD+ Tray Material Fast / Elite HD+ Maxi Tray Material Fast / Elite HD+ Monophase Normal / Elite HD+ Regular Body Normal / Elite HD+ Light Body Fast / Elite HD+ Light Body Normal / Elite HD+ Super Light Body Fast: vinyl polysiloxane, polydimethylsiloxane, silica fillers, methylhydrogenosiloxane, surfactant, pigments, inorganic aluminium fillers, organoplatinum complex.

5. COMPATIBLE DEVICES

- Zhermack 5:1 cartridges are compatible with the most common mixing machines available on the market, e.g. Sympress, Pentamix 2 and 3. Cartridges are not compatible with Pentamix Lite. Check compatibility with the manufacturer's instructions.
- The material is compatible with a range of dental stones for study model and working models (type 3 gypsum e.g. Zhermack's Elite Model or type 4 e.g. Zhermack's Elite Rock).

6. INDICATIONS FOR USE

The devices are intended to be used by dental professionals (user licensed to operate on the patient) for taking dental impressions on patients in a professional dental environment (public or private health facilities with the necessary regulatory approvals).

7. CONTRAINDICATIONS

Do not use on patients with known hypersensitivity to any of the components.

8. SIDE EFFECTS

Irritation, redness or signs of hypersensitivity may occur in case of allergy to any of the components.

9. CLINICAL BENEFITS

Impression material to accurately record the dimensions of tissues and/or oral prosthetic components and their spatial relationships.

10. PRECAUTIONS

- Insufficient data exist to support the use of VPS for impressions on exposed bone. VPS and accessory tips should only be used on or around intact mucosa.
- Devices marked "single use" on the labelling are intended for single use only. Do not reuse in other patients in order to prevent cross-contamination.
- Use of latex or "sulfur-containing" polymeric gloves may interfere with setting. Do not touch cords and prepared tooth/transfer with gloves. Rinse preparation with plenty of water before making the impression.
- Use of some hand soaps and lotions may interfere with setting.
- Never interchange different lots of base and catalyst, dosing spoons and jar lids to guarantee the working and setting times indicated in Table 14.
- Do not reuse the original cartridge cap. Store the cartridge with the disinfected used mixing tip.
- All products should be used at room temperature (23°C/73°F) to guarantee the working and setting times indicated in Table 14; higher temperatures reduce these times, lower temperatures increase them.
- To obtain optimum physical properties, do not use if the material is not homogeneously mixed.
- Check impression in the mouth (not on the bench) to be sure it is resistant before removal.
- Fitted and rigid trays with retentive features, pretreated with tray adhesive, are recommended. To reduce gagging, choking or swallowing, do not overfill tray.
- Material should dispense easily: DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated dispensing of the material. If the material doesn't dispense easily, discard the mixing tip, level component (bleed cartridge) and install a new tip.
- Do not use the material if the batch number and expiry date are missing or illegible and/or if package is broken. Do not use beyond the expiry date.
- If used with a syringe refer to the manufacturer instructions to check compatibility of the syringe with VPS before. Do not use a syringe that has been used with another chemical category.
- The colour of the tip must be the same as the colour of the cartridge cap to avoid excessive dispensing force, slower dispensing speed, incomplete mix and work/set time irregularities.
- Two-phase technique is NOT recommended for double arch trays.
- When used as a primary impression, material should be covered with a polyethylene spacer sheet, and protected from direct contact with oral tissue.
- Check the integrity of the impression and ensure that there is no residue of material to be removed from the patient's mouth. If necessary, rinse the mouth after impression removal.
- The two components of the device (base and catalyst) must be used in equal parts to guarantee the working and setting times indicated in Table 14. A higher dose of catalyst does not accelerate setting times.
- Use clean and properly disinfected and sterilised, or ready to use accessories.
- Replace the tip lock if it becomes visibly worn or damaged.
- Astringents based on aluminium or iron salts, hydrogen peroxide-containing dental materials, acrylic restoration resins and their residues may interfere with the setting reaction of material. If so, remove any residue before making the impression.
- Material should not be used with non-original accessories, intermixed or used in conjunction with other manufactur-

er's vinyl polysiloxanes or any polyether, polysulfide or conventional (condensation cured) silicones.

11. WARNINGS

- Causes eye irritation. Avoid direct contact of the product with the eyes. Protect eyes with airtight glasses. Protective eyewear is recommended for patients, too. In case of accidental contact with eyes, rinse immediately for several minutes with plenty of water and if irritation persists, seek medical attention.
- Avoid contact of individual components (base or catalyst) with skin to prevent irritation and possible allergic reaction. If contact with skin occurs immediately remove material with cotton and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- Avoid ingestion/swallowing of material to prevent irritation or obstruction. If accidental swallowing occurs, drink plenty of water. Seek medical attention in the event of digestive disorders.
- Avoid use on teeth with periodontal mobility and/or significant undercuts to prevent tooth extraction, injury to dental tissue or decementation of prosthetic artefacts.
- The product is intended for transient contact. Do not use in applications other than those for the intended use (e.g. temporary reliner) because it may cause irritation.

12. STEP-BY-STEP INSTRUCTIONS

12.1 PRELIMINARY OPERATIONS

Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients, too.

12.2 PREPARING THE DEVICES

Handle and dose outside the operating area and well away from splashes, sprays and aerosols of body fluids:

A. 250 / 450 ml JARS (Elite HD+ Putty Soft)

1. Remove seals and lids (including liner discs and films where present)
2. Take base and catalyst in equal parts. For a correct use it is recommended to use the white spoon for the catalyst and the coloured spoon for the base (fig. 3.1).
3. Close the jars. Do not swap over the lids and liner discs (where present).
4. Mix together equal amounts of Base and Catalyst, **until you obtain an evenly coloured mass without streaks** (fig. 3.2). Comply with the working time indicated in the technical data table (fig. 3.3).
5. Proceed to use the device on the patient (see paragraph 12.3).

B. 90 ml TUBES (Elite HD+ Regular Body / Elite HD+ Light Body)

1. Remove the base cap. Dispense a strip of base onto a clean surface. Close the base tube immediately after dispensing the material.
2. Open the catalyst tube.
3. Dispense a strip of catalyst the same length as the strip of base onto the clean surface (fig. 4.1). Repeat the procedure for each strip of base dispensed.
4. Close the catalyst tube immediately.
5. Mix together equal amounts of Base and Catalyst using a spatula for silicones, until you obtain an evenly coloured mass without streaks (fig. 4.2). Comply with the working time indicated in the technical data table.
6. Using a spatula and/or a syringe for use with elastomers, apply a small amount of the mixture where required, in accordance with the chosen technique (fig. 4.3 and 4.4).
7. Clean, disinfect and, if applicable, sterilise the spatula according to the instructions provided by its manufacturer.
8. Proceed in accordance with the chosen technique (see paragraph 12.3).

C. 1:1 50 ml CARTRIDGES (Elite HD+ Tray Material / Elite HD+ Monophase / Elite HD+ Regular Body / Elite HD+ Light Body / Elite HD+ Super Light Body)

1. Push the black lever underneath the rack upwards and pull the dispenser rack as far as it will go (fig. 5.1).
2. Lift the cartridge lock (fig. 5.1).
3. Insert the cartridge making sure that it fits the right direction. The cartridge must be inserted with the "V" on its base facing downwards (fig. 5.2). The cartridge lock can now be lowered easily.
4. After placing the cartridge on the dispenser, turn the cap by a 1/4 turn, pull it out with a firm movement and discard it (fig. 5.3).
5. Press the dispenser lever gently to dispense a small amount of the two components until they both flow out at the same time (fig. 5.4). **Remove any excess, taking care not to contaminate the two outlet holes.**
6. Select the appropriate mixing tip; the colour of the tip must be the same as the colour of the cartridge cap.
7. Insert the mixing tip in a clean motion so that the V-shaped notch on the coloured ring aligns with its housing on the cartridge (fig. 5.5). Turn the tip clockwise by a 1/4 turn, exerting pressure only on the coloured ring and not on the mixing tip itself (fig. 5.6).
8. Press the dispenser lever gently to dispense an appropriate amount of material. Do not use if the mix is not homogeneous. In this case, remove the dispenser from the operating field and remove the mixing tip; then dispense a small amount of the material to level the component and insert a new mixing tip.
9. Before use on each new patient, always apply a disposable cover/barrier system to the dispenser and cartridge, to be removed after each use (fig. 5.7). The disposable cover/barrier system does not replace the recommended cleaning and disinfection operations.
10. Proceed to use the device on the patient (see paragraph 12.3).

D. 5:1 380 ml CARTRIDGES (Elite HD+ Maxi Putty Soft; Elite HD+ Maxi Tray Material)

1. Press in the levers of the cartridge cap, remove the lid and discard it (fig. 6.1).

2. Switch on the 5:1 mixing machine and insert the cartridge, following the machine manufacturer's instructions.
3. Dispense a small amount of the two components until they both flow out at the same time. Remove any excess, taking care not to contaminate the two ports (fig. 6.2).
4. Insert the mixing tip, aligning the holes in the tip with the holes in the cartridge and fitting the hexagonal metal matrix of the mixing machine into the hexagonal matrix of the mixing tip (fig. 6.3).
5. Insert the tip lock, securing its fastenings correctly to the cartridge (fig. 6.4).
6. It is recommended to wrap the automatic mixer in disposable cellophane, removing it after each use.
7. Start the mixing machine and wait for the material to flow out. Always discard the first part of the material (fig. 6.5) and proceed according to the chosen technique, filling the impression tray (fig. 6.6) and/or elastomer syringe (fig. 6.7).
8. Proceed to use the device on the patient (see paragraph 12.3).

12.3 IMPRESSION TAKING PROCEDURE

Two-phase technique (double impression):

1. Apply a specific adhesive for addition silicones to the impression tray (following manufacturer's instructions).
2. Place a suitable amount of high viscosity material on the impression tray (fig. 1.1).
3. Place the loaded impression tray inside the patient's mouth within the working time (see paragraph 14) (fig. 1.2).
4. Remove the first impression from the patient's mouth when it is firmly hardened (for setting time see paragraph 14).
5. Wash the first impression and dry it thoroughly; then process it.
6. Apply the low-viscosity material where necessary (preparations, first impression, etc.) (fig. 1.3 and 1.4) and place the impression tray back inside the patient's mouth to take the second impression within the working time (see paragraph 14) (fig. 1.5).
7. Remove the impression from the patient's mouth when it is firmly hardened (for setting time see paragraph 14).
8. Clean and disinfect the impression before shipping and/or casting (see paragraph 12.4) and the dispenser, cartridge and mixing tip (see paragraph and 12.5).

One-step technique (simultaneous):

1. Apply a specific adhesive for addition silicones to the impression tray (following manufacturer's instructions).
2. Place simultaneously the higher viscosity material on the impression tray (fig. 2.1) and the lower viscosity material where necessary in the patient's mouth (preparations, transfer, etc.) (fig. 2.2).
3. Place the impression tray inside the patient's mouth within the working time (see paragraph 14) (fig. 2.3).
4. Remove the impression from the patient's mouth when it is firmly hardened (for setting time see paragraph 14).
5. Clean and disinfect the impression before shipping and/or casting (see paragraph 12.4) the dispenser, cartridge and mixing tip (see paragraph and 12.5).

12.4 MANUAL CLEANING AND DISINFECTION OF THE IMPRESSION

PRECAUTIONS / MANUAL CLEANING AND DISINFECTION OF THE IMPRESSION:

- Manual cleaning alone is not sufficient to condition an impression. Cleaning must always be followed by disinfection.
- Impressions must always be disinfected manually using a specific immersion disinfectant for addition silicones (vinyl polysiloxanes), (concentrated disinfectant based on quaternary ammonium salts and synergic co-formulants) or spraying a ready for use disinfectant, based on alcohol and surface tension reducers. Follow the chosen disinfectant manufacturer's instructions for contact time.
- The properly cleaned and disinfected impression must be stored in a closed, dry, clean environment well away from splashes, sprays and aerosols of body fluids and at an ambient temperature of around 23°C / 73°F.

12.4.1 PREPARATION PRIOR TO CLEANING AND DISINFECTION

1. Wear a pair of disposable gloves, a face mask, goggles and work clothes.
2. For disinfecting the impression by immersion (see paragraph 12.4.3 A), prepare a specific cleaning and disinfectant solution for impressions, approved according to local regulations and based on quaternary ammonium compounds and synergic co-formulants (e.g. Zeta 7 Solution), in a basin.

12.4.2 MANUAL CLEANING

1. Within a maximum of 5 minutes after removing the impression from the mouth rinse it under running water for **30 seconds** (fig. 7.1). If necessary, increase rinsing time until no dirt remains visible.
2. Disinfect the impression by immersion (see paragraph 12.4.3 A) or spray (see paragraph 12.4.3 B).

12.4.3 MANUAL DISINFECTION OF IMPRESSIONS

A. IMMERSION (fig. 7.2 A)

1. After cleaning (see paragraph 12.4.2), immerse the impression in the cleaning-disinfectant solution (see step 2 of paragraph 12.4.1) for the specified contact time.
2. Rinse and/or dry as instructed.

B. SPRAY (fig. 7.2 B)

1. After cleaning (see paragraph 12.4.2), place the impression in a transparent bag.
2. Insert the nozzle of the spray bottle of cleaner-disinfectant for impressions into the bag. Make sure that the spray is approved according to local regulations and based on ethanol, isopropyl alcohol and co-formulants (e.g. Zeta 7 Spray).
3. Holding the bag closed, spray the cleaner-disinfectant on all impression and impression tray's surfaces, and respecting the contact times.
4. Keep the bag closed and leave the product to act.
5. Remove the impression from the bag and dispose of the bag.

6. Rinse and/or dry as instructed.

12.5 MANUAL CLEANING AND DISINFECTION OF THE CARTRIDGE AND MIXING TIP

WARNINGS / PRECAUTIONS CONCERNING MANUAL CLEANING AND DISINFECTION OF THE CARTRIDGE AND MIXING TIP:

- Protective barriers are designed for single use and must be disposed of after each use in accordance with local regulations.
- Used cartridge and mixing tips must be processed as described below before storage. Manual cleaning alone is not sufficient for proper handling of the cartridge with the mixing tip inserted. Keep used, processed mixing tip on cartridge as a cap.
- Use only a disinfectant solution which is approved for its efficacy and use in accordance with the IFU of the disinfectant solution manufacturer.
- Perform manual cleaning and disinfection of the cartridge with the mixing tip attached within 1 hour of use.
- Before use, remove the used mixing tip from the cartridge and discard it. Level the components in the cartridge and fit a new mixing tip.
- The end of the life cycle of the cartridges is normally determined by the expiration and/or by the wear of the material herein.

12.5.1 DISASSEMBLY OF THE DEVICE AFTER USE

1. In the area used for cleaning and disinfection, remove the oral tips from mixing tips and the impression material residues inside.
2. Remove the disposable cover/barrier system.
3. Remove the contaminated gloves, sanitise hands and put on a clean pair of gloves.
4. Remove the cartridge with the mixing tip attached from the dispenser, taking care to avoid direct contact with the contaminated tip.
5. Clean and disinfect the dispenser according to the manufacturer's instructions.

12.5.2 MANUAL CLEANING

1. Do not remove the mixing tip from the cartridge.
2. Clean using disposable wipes completely impregnated with a disinfecting and cleaning solution*.
 - a. Take a disposable wipe completely impregnated with a disinfecting and cleaning solution* and wipe it **twice over all surfaces of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
 - b. Take a new disposable wipe completely impregnated with a disinfecting and cleaning solution* and wipe it **twice over all surfaces of the cartridge mixing tip**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
- Use additional wipes if necessary.
3. Use a clean, soft brush (e.g. Interlock® REF 09084) wetted with a disinfecting and cleaning solution* to clean the parts of the cartridge and the mixing tip with recesses and/or grooves (e.g. connection between the mixing tip and the head of the cartridge) taking care to remove all visible residues. Do not scrub the cartridge label with the brush.
4. Remove the disinfecting and cleaning solution* residues with a wet cloth soaked in water.

Proceed with the disinfection.

12.5.3 MANUAL DISINFECTION

1. Do not remove the mixing tip from the cartridge.
2. Disinfect using disposable wipes completely impregnated with a disinfecting and cleaning solution*.
 - a. Take a disposable wipe completely impregnated with a disinfecting and cleaning solution* and wipe it **twice over all surfaces of the cartridge**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
 - b. Take a new disposable wipe completely impregnated with a disinfecting and cleaning solution* and wipe it **twice over all surfaces of the cartridge mixing tip**, taking care to remove all visible residues and making sure that the fluid penetrates into all the recesses.
- Use additional wipes if necessary.
3. Allow the device to dry naturally for 5 minutes.
4. Inspect the cartridge and tip to make sure that all visible contamination has been removed. If visible contamination is still present, repeat the cleaning and disinfection steps described in points 12.5.2 and 12.5.3.
5. Before reusing the cartridge, insert it into the dispenser again, remove and dispose the used mixing tip (see paragraph 12.2).

The instructions provided above have been validated by the manufacturer of the medical device as the supplier of a reusable medical device. However, the responsibility for the cleaning and disinfection of the medical device lies with the person who performs the process, who must guarantee that it is carried out correctly in the facility where cleaning and disinfection take place and that the desired result is achieved. This requires verification and/or validation and continuous monitoring of the process.

*bactericidal, yeasticidal, tuberculocidal and virucidal disinfecting and cleaning solution containing mixtures of ethanol (≈35%) and isopropyl alcohol (≈15%) (e.g. Zeta 3 Soft, disinfecting and cleaning solution, with a contact time of 5 minutes) approved in accordance with local regulations; use it in compliance with the instructions provided by the manufacturer of the disinfectant solution.

13. INFORMATION ON THE CORRECT USE OF PACKAGING AND RELATED ACCESSORIES

Use accessories that are clean and properly disinfected and sterilised or ready to use.

The packaging of the products and related accessories, excluding those mentioned in section 12.5, can be re-used safely without special cleaning and disinfection processes provided the following recommendations for preservation, dosing and storage are observed:

- Handle, dose from and store packaging and related accessories outside the operating area, in clean facilities and well away from splashes, sprays and aerosols of body fluids;
- Always wear new, uncontaminated gloves when handling jars and dosing spoons. If gloves become contaminated, change them immediately before handling the product;
- Close packaging immediately after dosing.

Dispose of contaminated products and related accessories as instructed in disposal paragraph 15.

14. TECHNICAL DATA

1 - ISO 4823

2 - Mixing ratio Base: Catalyst

3 - Mixing time

4 - Working time (including mixing time)*

5 - Time in the oral cavity (min.)

6 - Setting time*

7 - Elastic recovery

8 - Linear dimensional change 24 h

9 - Shore-A hardness 1 hour

*The times mentioned are intended from the start of mixing at 23°C / 73°F.

15. STORAGE, STABILITY AND DISPOSAL

The devices should be used within a professional dental environment (legally-certified public or private health facilities).

- Dispose of cartridge when empty. Do not attempt to refill used cartridges. Oral tips with material inside must be disposed of as special waste at risk of biological contamination.
- Based on the information in the SDS, the waste generated by the disposal of the product does not represent physical danger, for the person or the environment. It can therefore be managed, according to local regulations, as a non-dangerous waste. For the manipulation of contaminated components refer to the paragraphs 12.4 and 12.5. If the device and the impression are contaminated, dispose of them as special waste characterised by a risk of biological contamination. If they are not contaminated, dispose of them according to applicable local regulations.

Do not use after the expiry date. Store the product at a temperature between 5°C / 41°F and 27°C / 80°F. Do not store the product in direct sunlight.

16. IMPRESSION CASTING

- The impression may be poured in plaster 60 minutes after disinfection or delayed up to 15 days.
- The material is compatible with a range of dental stones for study model and working models (type 3 gypsum e.g. Zhermack's Elite Model or type 4 e.g. Zhermack's Elite Rock).

17. CONSERVATION OF THE IMPRESSION

Store the impressions at ambient temperature, in a dry place away from direct sunlight. If the impression is to be shipped, use suitable packaging to prevent distortion.

18. IMPORTANT REMARKS

Information provided in any way, even during demonstrations, does not invalidate the instructions for use. Operators are required to check that the product is suitable for the application envisaged. The manufacturer cannot be held responsible for damage, including to third parties, deriving from failure to follow instructions or from unsuitability for an application. The manufacturer's liability is, in any case, limited to the value of the products supplied. Report any serious incident involving the medical device to the manufacturer and to the relevant authorities.

Three-year shelf-life from the production date.

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3. VERPAKKING

- Potten 250 ml basis + 250 ml katalysator (Elite HD+ Putty Soft)
- Potten 450 ml basis + 450 ml katalysator (Elite HD+ Putty Soft)
- Patroon 380 ml (Elite HD+ Maxi Putty Soft; Elite HD+ Maxi Tray Material)
- Patroon 50 ml (Elite HD+ Tray Material; Elite HD+ Monophase; Elite HD+ Regular Body; Elite HD+ - Light Body; Elite HD+ Super Light Body)
- Tube 90 ml (Elite HD+ Regular Body; Elite HD+ Light Body)

4. SAMENSTELLING

Elite HD+ Putty Soft Normal / Elite HD+ Maxi Putty Normal / Elite HD+ Putty Soft Fast /

Elite HD+ Maxi Putty Fast: vinylpolysiloxaan, polydimethylsiloxaan, siliciumdioxidevulstoffen, methylwaterstofsiloxaan, koolwaterstoffen, pigmenten, anorganische aluminiumvulstoffen, additieven, organoplatinacomplex.

Elite HD+ Tray Material Fast / Elite HD+ Maxi Tray Material Fast / Elite HD+ Monophase Normal / Elite HD+ Regular Body Normal / Elite HD+ Light Body Fast / Elite HD+ Light Body Normal / Elite HD+ Super Light Body Fast: vinylpolysiloxaan, polydimethylsiloxaan, siliciumdioxidevulstoffen, methylwaterstofsiloxaan, oppervlakteactieve stof, pigmenten, anorganische aluminiumvulstoffen, organoplatinacomplex.

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- 450 krukke med base + 450 ml katalysator (Elite HD+ Putty Soft)

- 380 ml magasin (Elite HD+ Maxi Putty Soft; Elite HD+ Maxi Tray Material)

- 50 ml magasin (Elite HD+ Tray Material; Elite HD+ Monophase; Elite HD+ Regular Body; Elite HD+ - Light Body; Elite HD+ Super Light Body)

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Elite HD+ Putty Soft Normal / Elite HD+ Maxi Putty Normal / Elite HD+ Putty Soft Fast / Elite HD+ Maxi Putty Fast: vinylpolysiloxan, polydimethylsiloxan, silicafyldstoffer, methylhydrogenosiloxan, kulbrinter, pigmenter, uorganiske aluminiumfyldstoffer, additiver, organoplatinum complex.

Elite HD+ Tray Material Fast / Elite HD+ Maxi Tray Material Fast / Elite HD+ Monophase Normal / Elite HD+ Regular Body Normal / Elite HD+ Light Body Fast / Elite HD+ Light Body Normal / Elite HD+ Super Light Body Fast: vinylpolysiloxan, polydimethylsiloxan, silicafyldstoffer, methylhydrogenosiloxan, overfladeaktivt stof, pigmenter, uorganiske aluminiumfyldstoffer, organoplatinum complex.

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TECHNICAL DATA										
	Elite HD+	Putty Soft Putty Soft Maxi Normal Set	Putty Soft Putty Soft Maxi Fast Set	Tray Material Fast Set	Tray Material Maxi Fast Set	Monophase Normal Set	Regular Body Normal Set	Light Body Normal Set	Light Body Fast Set	Super Light Body Fast Set
1	ISO 4823	Type 0 Putty Consistency	Type 0 Putty Consistency	N/A	Type 1 Heavy-bodied Consistency	Type 2 Medium-bodied Consistency	Type 2 Medium-bodied Consistency	Type 3 Light-bodied Consistency	Type 3 Light-bodied Consistency	Type 3 Light-bodied Consistency
2	Mixing ratio Base: Catalyst	1:1 5:1	1:1 5:1	1:1	5:1	1:1	1:1	1:1	1:1	1:1
3	Mixing time	In jar 30" Auto mixing System	In jar 30" Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	In tube 30" Auto mixing System	In tube 30" Auto mixing System	Auto mixing System	Auto mixing System
4	Working time (including mixing time)*	2'	1'30"	1'30"	1'30"	2'	2'	2'	1'30"	1'30"
5	Time in the oral cavity (min.)	3'30"	2'30"	2'30"	2'30"	3'30"	3'30"	3'30"	2'30"	2'30"
6	Setting time*	5'30"	4'	4'	4'	5'30"	5'30"	5'30"	4'	4'
7	Elastic recovery	99%	99%	99%	99%	99%	99%	99%	99%	99%
8	Linear dimensional change 24 h	0,3%	0,3%	0,3%	0,3%	0,3%	0,3%	0,3%	0,3%	0,3%
9	Shore-A hardness 1 hour	65	65	70	65	65	55	55	55	50