Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

PRODUCTS COVERED UNDER THIS IFU:

OverdenSURE® Abutments and components.

ELABELING:

The Instructions for Use can be accessed online by visiting the website provided on the label. Additional translations are also available in electronic format for download.

CAUTION – PLEASE READ CAREFULLY

Products distributed by ZimVie US Corp LLC should only be used by dental specialists with experience in maxillofacial implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. Products marked for single use exhibit a "Do not re-use" symbol on the product label. Products bearing this symbol on the label must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health such as patient tissue infection.

Products covered in this document must be dry fitted before use to check that the fit is correct. The clinician is responsible for the correct application of restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor.

During any intraoral use and manipulation all products must be secured to prevent aspiration due to its shape and size.

Caution: U.S. Federal law restricts this device to be sold to, or on the order of, a licensed dentist or physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must consist of a screw check. If the screws are subject to unusual wear, the complete integrity of the implant abutment should be checked. Failure to follow this instruction puts the patient at risk and will void the warranty.

INTENDED USE AND INSTRUCTIONS FOR USE - PLEASE READ CAREFULLY

The OverdenSURE® Attachment System is designed for use with overdentures or partial dentures, retained in whole or in part by endosseous implants in the mandible or maxilla.

OverdenSURE Attachment Placement Indications

IMPLANT SYSTEM	PLATFORM	INTRAORAL POSITION	INDICATION
Eztetic™ 3.1mm(D) TSX™ 3.1mm(D)	NP - 2.9mm(D)	Anterior	Yes
		Pre-molar	Splinted
		Molar	Not indicated
Tapered Screw-Vent® Trabecular Metal™ TSX™	3.5mm(D) 4.5mm(D)	All	Yes
Tapered Screw-Vent® Trabecular Metal™	5.7mm(D)	All	Yes
Certain® External Hex	3.4mm(D) 4.1mm(D) 5.0mm(D)	All	Yes
Spline®	3.25mm(D) 3.75mm(D) 5.0mm(D)	All	Yes
Tapered SwissPlus®	3.8mm(D) 4.8mm(D)	All	Yes
ANKYLOS® C/X (DA-INTERNAL)	3.5/4.5mm(D)	All	Yes
Astra Tech Implant System® Osseospeed® EV (DAEV-CONICAL)	3.6mm(D) 4.2mm(D) 4.8mm(D)	All	Yes
Astra Tech Implant System® Osseospeed® TX (DATX-INTERNAL)	3mm(D) 3.5/4mm(D) 4.5/5mm(D)	All	Yes
Biohorizons Tapered Internal (BHI-INTERNAL)	3.5mm(D) 4.5mm(D)	All	Yes
Bone Level (SBL-CONICAL)	3.3mm(D) 4.1/4.8mm(D)	All	Yes
CAMLOG® (CCA-INTERNAL)	3.8mm(D) 4.3mm(D)	All	Yes
MIS® SEVEN® Internal Hexagon (MSH-INTERNAL)	3.75/4.2mm(D) 5/6mm(D)	All	Yes
Nobel Brånemark System® (NB-EXTERNAL)	3.5mm(D) 4.1mm(D)	All	Yes
Nobelactive® & Nobelreplace® Conical (NAC-CONICAL)	3.5mm(D) 4.3/5mm(D)	All	Yes
Nobelreplace® & Replace Select™ Tapered (NRT-TRILOBE)	3.5mm(D) 4.3mm(D)	All	Yes
TISSUE LEVEL (STL-TISSUE)	4.1/4.8mm(D)	All	Yes
XIVE® (DX-INTERNAL)	3.4mm(D) 3.8mm(D)	All	Yes
Neodent GM (NG-CONICAL)	3.0mm(D)	All	Yes

CONTRAINDICATIONS – PLEASE READ CAREFULLY

All materials used are biocompatible; however, some patients may be allergic or hypersensitive to any of the materials and their components

- Implant Attachment: Not appropriate where a totally rigid connection is required.
- Use of a single implant with divergence of greater than 20 degrees from vertical is not recommended.

The 2.9mm platform OverdenSURE Attachments should not be used in the molar region

INTENDED USERS AND PATIENT GROUPS

Products covered in this document should only be used by dental specialists with experience in dental implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. Their use is limited to dental laboratories and clinics.

Indicated for edentulous patients (totally or partially) who require oral rehabilitation by means of implant-supported prostheses. Rehabilitation can be unitary, multiple or overdentures, both in the upper and lower jaw. Its use is indicated in patients with complete dentofacial development. There are no differences in end users by age group, sex, ethnicity, family predisposition or genetic aspects. It is not indicated for use in patients without dental problems.

CLINICAL BENEFITS AND UNDESIRABLE SIDE EFFECTS

As a clinical benefit, patients can expect full or partial rehabilitation of the dentition, allowing them to regain proper masticatory function.

No side effects directly related to the use of ZimVie implantology products have been described, however, the use of these devices may be part of an invasive treatment that may be associated with typical side effects such as inflammation, bleeding, hematoma, pain or swelling.

SERIOUS INCIDENT NOTICE

For patients / users / third parties in the European Union with an identical regulatory regime (Regulation 2017/745/EU) if, due to the use of the product, a serious incident occurs, notify ZimVie US Corp LLC at:

ZimVie US Corp LLC – Headquarters
Address: 4555 Riverside Drive, Palm Beach Gardens, Florida 33410
E-mail: DentalInternationalComplaints@ZimVie.com

WARNINGS – PLEASE READ CAREFULLY

Products distributed by ZimVie US Corp LLC have not been evaluated for safety and compatibility in the Magnetic Resonance environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of products distributed by ZimVie US Corp LLC in the Magnetic Resonance environment is unknown. Scanning a patient who has this device may result in patient injury.

Magnetic Resonance Imaging (MRI) Safety Information

MR Conditional

Warning: The RF safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.

A patient with this device can be scanned safely in an MR system under the following conditions:

Static Magnetic Field Strength (B0)	1.5 T and 3.0 T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil. Extremity T/R coils permitted. Excludes Head T/R coil.
Operating Mode	Normal Operating Mode in the allowed imaging zone
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head specific absorption rate (SAR)	Not evaluated for head landmark
Scan Duration	No specific constraints due to implant heating

Note: the removable restorations should be taken out prior to scanning.

STERILIZATION AND RE-USE

Products covered in this document are supplied NON-STERILE. See table 1. The recommended sterilization cycle is standard gravity autoclave, exposure at 121 °C / 250 °F for 30 minutes with a drying time of 30 minutes (in accordance with the UNE-EN ISO 17665-1 and -2), using a sterilization wrap that is EU and FDA cleared for the indicated cycle.

Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Discard single-use and worn reusable instruments immediately after use. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

Frequent use has a minimal effect on the instruments. The end of instrument life is usually determined by wear and tear and deterioration during use. Therefore, instruments can be reused with due care, provided they are not damaged or contaminated.

For the reuse of some instruments, they must be cleaned prior to sterilization. Here are some cleaning tips:

- Never place instruments of different types of materials together.
- Do not use metal brushes to remove impurities.
- Use disposable syringes for instrument cavities.
- When selecting detergents and disinfectants, make sure they are commercially available products and always follow the manufacturer's instructions.

SHELF LIFE

The abutments are intended to be handled and manipulated by professional dental technicians before being shaped into their final form. They are labelled "for single use", but do not degrade over time and are not susceptible to variations in environmental or storage conditions. Therefore, neither the shelf life of the device nor the shelf life of the package is limited and no "use by" date is provided on the label.

Clinical screws are intended to be handled by professionals for implant-prosthesis fixation, and once assembled, they are not susceptible to variations in environmental conditions. Their shelf life is not limited.

DISPOSAL

Disposal of devices should follow local regulations and environmental requirements, taking into account different levels of contamination.

Table 1 General Information

FAMILY	COMMERCIAL CODES		RAW MATERIAL	STERILIZATION	RE-USE
OverdenSURE® Abutment/ Class IIb	ODS-CERTXXX ODS-EXHEXXX ODS-TSVXXX ODS-NPXXX ODS-SPLXXX ODS-SWPXXX ODS-DAXXX ODS-DAEVXXX ODS-DATXXX ODS-BHIXXX ODS-BHISXXX ODS-SBLXXX ODS-CCAXXX ODS-MSHXXX ODS-NBXXX ODS-NACXXX ODS-NRTXXX ODS-STLXXX ODS-DXXX ODS-NGXXX	ODSA-CERTXXX, ODSA-EXHEXXX ODSA-TSVXXX ODSA-NPXXX ODSA-DAXXX ODSA-DAEVXXX ODSA-DATXXX ODSA-BHIXXX ODSA-BHIAXXX ODSA-SBLXXX ODSA-CCAXXX ODSA-MSHXXX ODSA-NBXXX ODSA-NACXXX ODSA-NRTXXX ODSA-STLXXX ODSA-DXXX ODSA-NGXXX (ODSA product codes available only in U.S.A)	Titanium alloy ELI Ti-6Al-4V, ZrN coated	Autoclave before intraoral use	Single use
OverdenSURE® Abutment Retention Components/ Class IIa	ODS-RIKIT2001 ODS-RIKIT2002 ODS-RIKIT4001 ODS-RIKIT4002 ODS-RIR4004 ODS-RIB2004 ODS-RIP2004 ODS-RIC2004 ODS-RIY4004 ODS-RIG4004	ODSZ-RIKIT2001 ODSZ-RIKIT2002 ODSZ-RIKIT4001 ODSZ-RIKIT4002 ODSZ-RIR4004 ODSZ-RIB2004 ODSZ-RIP2004 ODSZ-RIC2004 ODSZ-RIY4004 ODSZ-RIG4004 (ODSZ product codes available only in U.S.A)	Retention Insert - Polynil (Nylon)	The nylon Retention Inserts may be sterilized/ disinfected using an FDA approved liquid chemical sterilant. Nylon Retention Inserts must be soaked for a minimum of 3 hours in the liquid sterilant at room temperature	Single use
OverdenSURE® Abutment Retention Components/ Class I	ODS-RH01 ODS-RH04 ODS-RH010 ODS-BS020		Retention Housing – Titanium Alloy (Ti-6Al-4V ELI) Blockout Spacer – White Silicone	N/A	Single use
OverdenSURE® Analog	ODS-AA01 ODS-AA04 ODS-AA010		Surgical AISI-303	No sterilization needed. Not indicated for intraoral use.	Not recommended
OverdenSURE® Impression Coping	ODS-IC01 ODS-IC04 ODS-IC010		Surgical AISI-303 Screw: Titanium Alloy Ti-6Al-4V ELI	Autoclave before intraoral use	Not recommended
OVERDENSURE® Bar Attachment	ODS-BARATT01 ODS-BARATT02 ODS-BARATT04		Titanium alloy ELI Ti-6Al-4V, ZrN coated	Autoclave before intraoral use	Single use
OverdenSURE® Multi-Unit	ARMU-ODA483		Titanium alloy ELI Ti-6Al-4V, ZrN	Autoclave before intraoral use	Single use

Abutment Attachment/ Class IIb		coated		
TOOLS / Class I	ODS-DRVR ODS-IRTOOL ODS-DRVRA ODS-ADX	Surgical stainless steel SS 316L Tips: stainless steel AISI 420 MOD	Autoclave before intraoral use	Reusable
OverdenSURE® Scan Body	ODS-SB010	Polyether ether ketone (PEEK)	Autoclave before intraoral use	Reusable

Table 2 Compatibility table.

OVERDENSURE

IMPLANT SYSTEM COMPATIBILITY	PLATFORM	GINGIVAL HEIGHT
Eztetic™ 3.1 mm(D) TSX™ 3.1 mm(D)	NP - 2.9mm(D)	1, 2, 3, 4, 5, 6mm(H)
Tapered Screw-Vent® Trabecular Metal™ TSX™	3.5mm(D) 4.5mm(D)	1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H)
Tapered Screw-Vent® Trabecular Metal™	5.7mm(D)	1, 2, 3, 4, 5, 6mm(H)
Certain® External Hex	3.4mm(D) 4.1mm(D) 5.0mm(D)	1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H)
Spline®	3.25mm(D) 3.75mm(D) 5.0mm(D)	1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H)
Tapered SwissPlus®	3.8mm(D) 4.8mm(D)	1, 2, 3, 4, 5, 6mm(H) 1, 2, 3, 4, 5, 6mm(H)
ANKYLOS® C/X (DA-INTERNAL)	3.5/4.5mm(D)	2, 3, 4mm(H)
Astra Tech Implant System® Osseospeed® EV (DAEV-CONICAL)	3.6mm(D) 4.2mm(D) 4.8mm(D)	1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H)
Astra Tech Implant System® Osseospeed® TX (DATX-INTERNAL)	3mm(D) 3.5/4mm(D) 4.5/5mm(D)	4mm(H) 1, 2, 3, 4mm(H) 1, 2, 3, 4, 5mm(H)
Biohorizons Tapered Internal (BHI- INTERNAL)	3.5mm(D) 4.5mm(D)	1, 2.5, 3.5, 4.5, 5.5mm(H) 1, 2.5, 3.5, 4.5mm(H)
Bone Level (SBL-CONICAL)	3.3mm(D) 4.1/4.8mm(D)	2, 3, 4mm(H) 1, 2, 3, 4mm(H)
CAMLOG® (CCA-INTERNAL)	3.8mm(D) 4.3mm(D)	1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H)
MIS® SEVEN® Internal Hexagon (MSH- INTERNAL)	3.75/4.2mm(D) 5/6mm(D)	1, 2.5, 3.5, 4.5, 5.5mm(H) 1, 2.5, 3.5, 4.5mm (H)
Nobel Brånemark System® (NB- EXTERNAL)	3.5mm(D) 4.1mm(D)	1, 2, 3mm(H) 1, 2, 3, 4, 5mm(H)
Nobelactive® & Nobelreplace® Conical (NAC-CONICAL)	3.5mm(D) 4.3/5mm(D)	1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H)
Nobelreplace® & Replace Select™ Tapered (NRT-TRILOBE)	3.5mm(D) 4.3mm(D)	1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H)
TISSUE LEVEL (STL-TISSUE)	4.1/4.8mm(D)	1, 2, 3, 4mm(H)
XIVE® (DX-INTERNAL)	3.4mm(D) 3.8mm(D)	1, 2, 3, 4mm(H) 1, 2, 3, 4mm(H)
Neodent GM (NG-CONICAL)	3.0mm(D)	0.8, 1.5, 2.5, 3.5, 4.5, 5.5mm(H)
Multi-Unit Abutment	4.8 mm	3mm(H)

INTENDED USE AND INSTRUCTIONS FOR USE - OVERDENSURE® ATTACHMENTS AND RETENTION COMPONENTS

OverdenSURE® Abutment is a resilient implant attachment used for prosthetic restorations prepared by dental technicians in dental laboratories or by dental professionals chairside. Recommended for removable dentures.

OverdenSURE® Abutment retention system is a part of the anchoring system for full or partial overdentures on implants, in cases where conventional full dentures have problems with retention and stability. The male retention kit can be mounted in the dental clinic using the self-curing resin, and in the dental laboratory, using the standard procedure.

Refer to the OverdenSURE® Attachment System Technique Manual for procedure instructions (ZBINST110).

Retention force/divergence between implants: blue (1.5 lbs/680g, up to 20°), pink (3 lbs/1360g, up to 20°), clear (5 lbs/2270g, up to 20°), black (Processing Nylon male, 0 lbs, 0g, up to 40°), red (1 lbs/450g, up to 40°), yellow (2 lbs/910g, up to 40°), green (4 lbs/1810g, up to 40°). ODS-RIKIT2001 and ODS-RIKIT2002 are kits up to 20°, ODS-RIKIT4001 and ODS-RIKIT4002 are kits up to 40°.

INTENDED USE AND INSTRUCTIONS FOR USE - TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. Exchangeable screwdriver tips may be used for various implant systems. For intraoral use, all tools must be thoroughly autoclaved before use. Where using removable heads, make sure that it cannot come loose from the screwdriver. The screwdriver should always be fastened with a safety pin and cord attached to the groove designed for this purpose. It should be handled with care, taking the necessary precautions not to hurt the patient. Before using the screwdriver, check that it is the correct size and shape for the screw head that will be used. Only use screwdrivers that are in perfect conditions. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

In case of the need for adaptors, please assure that they correspond to 4 x 4 connection and to its abutment system (OVERDENSURE®).

SUMMARY ON SAFETY AND CLINICAL PERFORMANCE (SSCP)

The summary of safety and clinical performance (SSCP) for products belonging to the families of devices is available in the European Database on Medical Devices (EUDAMED). You can consult the document linked to the basic UDI-DI (8435457210NS02F4N, 8435457212NS02A56) on EUDAMED's public website: <https://ec.europa.eu/tools/eudamed/#/screen/home>

SYMBOLS GLOSSARY

ANSI/AAMI/ ISO 15223-1:2012 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements.

Symbol	Title of Symbol	
	Consult Instructions for Use or Consult Electronic Instructions for Use labeling.zimvie.com	Indicates the need for the user to consult the instructions for use.
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Unique device identifier	
	Model number	

	Manufacturer	Indicates the medical device manufacturer.
	Date of Manufacture	Indicates the date when the medical device was manufactured
	Distributor	
	NON Sterile product	
	MR Conditional	Magnetic Resonance (MR) Safety Information
	Medical Device	
	European Conformity	
	European Conformity (Class IIb and Class IIa)	
Rx Only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	

All products distributed by ZimVie US Corp LLC are manufactured in accordance with ISO 13485 and MDSAP requirements and bear the CE marking.

These instructions replace all previous editions.

For any information about Products distributed by ZimVie US Corp LLC, please contact your local Territory Manager.



ZimVie US Corp LLC
4555 Riverside Drive
Palm Beach Gardens, FL
33410, USA
+1-561-776-6700
<http://www.zimvie.com/en/dental.html>

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ASTRA TECH™, OsseoSpeed®, ANKYLOS® C/X, and XiVE® are registered trademarks of DENTSPLY Implants.
BioHorizons® is a registered trademark of BioHorizons, INC.
Brånemark System®, NobelReplace® and Replace Select™, NobelActive®, NobelActive™, NobelReplace® CC are registered trademarks of Nobel Biocare AB.
MIS® and SEVEN® are registered trademarks of MIS Implant Technologies Ltd.
OverdenSURE®, Tapered Screw-Vent®, Eztetic™, TSX™, Trabecular Metal™, Certain®, Spline® and Tapered Tapered SwissPlus® are registered trademarks of ZimVie US Corp LLC.
Camlog® is a trademark of Camlog Biotechnologies Group.

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Brånemark System®, NobelReplace® ja Replace Select™, NobelActive®, NobelActive System™, NobelReplace® CC are registered trademarks of Nobel Biocare AB:n rekisteröityjä tavaramerkkejä.
For informasjon om produkter som distribueres av ZimVie US Corp LLC, vennligst kontakt din lokale områdeansvarlige.
MedNet SWISS GmbH. D4 Platz 4. 6039 Root D4. Sveits
MIS® og SEVEN® er registrerte varemerker som tilhører MIS Implant Technologies Ltd.