



STA Single Tooth Anesthesia[®] System New Trolley



Item Name: STA Single Tooth Anesthesia[®] System New Trolley

Order #: WA-8012

Shipping Quantity: 1 trolley requiring assembly

Dimensions & Weight: 25lbs, 30" x 17" x 7"

11.34 kg, 76.2 cm x 43.18 cm x 17.78 cm



Description:

- New STA trolley has improved rolling capability with wheels that have a high load capacity - ideal for use on carpeted and tile floors.
- Trolley uses lightweight steel with quality welds - assuring the Trolley will be durable and functional for many years to come.
- Top shelf has a screw mount that attaches to the bottom of the STA instrument - holding it securely on the Trolley and to keep it from falling or breaking.
- Trolley has a tray that will accommodate a standard box of STA handpieces - making accessing the disposables to use with each patient very easy.
- Handle at front to help pull and control the direction of the Trolley - making sure it arrives at the desired location with very little effort.

***** STA Instrument not included. The Trolley works with all Wand[®] Systems.**