

SCREW REMOVER KIT



1/ SR kit

- SR kit is consisted of various tools designed to remove the fracture screw quickly and safely. (1,200 RPM without irrigation)

2/ Key Points of Neo Screw Remover kit

- Remove fractured screw without damaging the fixture or pulling out the fixture.
- Fractured screw can be removed in 3 minutes.
- Various implant system can be applied to Neo Screw Remover kit.
- Screw with excessive torque can be removed with ratchet.

Screw Remover Kit Technique Guide




1 Fractured screw

2 Perfect Guide



3 Reverse Drilling

4 Screw Removing



Screw Remover Kit Clinical Case




Ordering information

Catalog Number	Description
SRKIT	Screw Remover kit

Note

SR kit consisting of the tools designed to remove the fracture screw fast and safely in case that the screw inside the inserted fixture are fractured due to several reasons



Neobiotech Co.,Ltd

104, 24 / 102, 103, 104-1, 104-2, 105, 106, 205, 212, 312, 509, 510, 511, 10F,
36 Digital-ro 27-gil, Guro-gu, Seoul, Korea Tel 82-2-582-2885, Fax 82-2-582-2883
www.neobiotech.com / Made in Korea

Dongbang Acuprime

1 Forrest Units, Hennock Road East, Marsh Barton, Exeter EX2 8RU, UK
Tel +44 1392-829500 Fax +44 1392-829232

**Neo
Biotech**
Satisfaction to Dentists