

Universal Prosthetic Kit

• The Universal Prosthetic Kit is composed of various prosthetic drivers that are compatible to different brands and sizes of implant abutments.

/Intended Use

• This Universal Prosthetic Kit is used when implanting a fixture in a place where a tooth is missing, or when connecting or detaching a cover screw or other upper structures.

Component

1. Driver

This tool is used to tighten an abutment or abutment screw.

1.2 Hex short, long / 1.25 Hex short, long / SCS short, long / Unigrip short, long driver

Name	Picture		Compatible implant or company	Code	Length(L)
1.2 Hex Driver (0.48'')		Neobiotech Biomet 3i	HD1210S	10mm	
			Dentsply-Friadent, Xive, Osstem, DIO, Megagen	HD1215S	15mm
1.25 Hex Driver (0.5")			Biohorizons Zimmer, Astra Implant Direct Legacy Camlog, MIS, Dentis, Dentium	HD12510S	10mm
				HD12515S	15mm
SCS Driver			Straumann Implant Direct Replant	SCSD10S	10mm
				SCSD15S	15mm
Unigrip Driver		-	Nobel Biocare Implant Direct Swishplait Southern Implants Neoss	UGD10S	10mm
				UGD15S	15mm

2. Driver Holder Code : DH01

The Driver Holder connects to various prosthetic drivers to allow easy handling of the healing cover screws and abutment screws.

3. Torque Ratchet Code : TW60

The Torque Ratchet is used to tighten the hex driver and measure torque. By changing the direction of the tip, tightening and loosening can be controlled.





